plusmed

health management

Model: **COMPACT/NEB**

COMPRESSOR NEBULIZER
KOMPRESÖRLÜ NEBULİZATÖR
NÉBULISEUR À COMPRESSEUR
COMPRESOR NEBULIZADOR

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EN Instruction Manual

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ES Manual de Instrucciones

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CONTENTS

IEC SYMBOLS	1
IMPORTANT SAFEGUARDS	1
INTRODUCTION	2
CONSTRUCTION	2
COMPRESSOR	2
REUSABLE NEBULIZER KIT	2
HOW TO OPERATE YOUR COMPRESSOR	3
NEBULIZER CLEANING	3
COMPRESSOR CLEANING	4
FILTER REPLACEMENT	4
MAINTENANCE.....	5
SPECIFICATIONS.....	5
GUIDELINE AND DECLARATION.....	6

IEC SYMBOLS

	Manufacturer		Date of manufacture
	Follow guidebook		Type-BF applied part
	EU Authorized representative		Class II Equipment
IP21	The compressor nebulizers are protected against solid foreign objects of >Φ12,5mm, and vertically falling water drops per IEC60529		Serial number
	The device contains electrical and/or electronic equipment that must be recycled per Directive 2012/19/EU – Waste Electrical and Electronic Equipment (WEEE)		

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed.

Read all instructions before using. Important information is highlighted by these terms:

DANGER-	Urgent safety information for hazards that will cause serious injury or death.
WARNING-	Important safety information for hazards that might cause serious injury.
CAUTION-	Information for preventing damage to the product.
NOTE-	Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

DANGER

To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place or drop into water or other liquid.
5. Do not reach for a product that has fallen onto water. Unplug immediately.

WARNING

To reduce the risk of burns, electrocution, fire or injury to persons:

1. The product should never be left unattended when plugged in.
2. Close supervision is necessary when this product is used by, on, or near children or physically incapacitated individuals.
3. Used this product only for its intended use as described in this guidebook. Do not use attachments not recommended by the manufacturer.
4. Never operate this product if:
 - a) It has a damaged power cord or plug.

- b) It is not working properly.
 - c) It has been dropped or damaged.
 - d) It has been dropped into water.
5. Return the product to an authorized service centre for examination and repair.
 6. Keep the power cord away from heated surfaces.
 7. Never block the air openings of the product or place it on a soft surface, such as a bed or couch, where the air openings may be blocked. Keep the air openings free of lint, hair, and the like.
 8. Never use while drowsy or asleep.
 9. Never drop or insert any object into any opening or hose.
 10. Do not use in an oxygen enriched environment.
 11. Do not operate unit in a backpack or carrying bag.

NOTE- *this unit employs a temperature and/or current limiting device that will shut off the entire product should a malfunction occur or in the event of product misuse, abuse, or abnormal operation.*

SAVE THESE INSTRUCTIONS

INTRODUCTION

Your doctor has prescribed a liquid medication to treat your respiratory disorder. To best utilize this liquid medication, he has prescribed a Plusmed CompactNeb compressor nebulizer. Your Plusmed CompactNeb compressor nebulizer works to convert the medication to a high-quality mist of fine particles that penetrates deep into the lungs.

Be sure that you read and understand the information in this guidebook. By following these simple instructions and the advice of your physician, the compressor will become an effective addition to your therapeutic routine.

Intended Use

The Plusmed CompactNeb compressor nebulizer includes an AC powered air compressor that provides a source of compressed air for both professional and home health care application. The compressor is used in conjunction with a jet nebulizer to convert certain inhaled drugs into an aerosol form for inhalation by patients with acute asthma. The device can be used for adult or paediatric patients. To reduce cross infection, the reusable nebulizer kit is intended for single patient use.

Your clinical response to the prescribed treatment shall be regularly supervised by your clinician. The most important determinant of appropriate dosing is the clinician's judgement of your response.

CONSTRUCTION

COMPRESSOR

REUSABLE NEBULIZER KIT

1. Mouth piece
2. Mask (adult)
3. Mask (paediatric)
4. Medication cup

5. Baffle
6. Air tube

HOW TO OPERATE YOUR COMPRESSOR

NOTE- before initial operation, your nebulizer should be cleaned following cleaning procedure, or as recommended by your doctor.

1. Make sure the power switch is in “Off”(O) position.
2. Place compressor so that the rubber feet are resting on a level, sturdy surface and upon which the controls can be easily reached when seated. Never block the air openings of the product or place it on a soft surface, such as a patient’s lap or a bed or couch where the air openings may be blocked.
3. Plug power cord into an appropriate wall outlet.

DANGER

Plusmed CompactNeb compressor nebulizer must be operated on the specified power source to avoid risk of electric shock and damage to the compressor.

4. Connect one end of the air tube to the compressor air-outlet port.

NOTE- during high-humidity weather conditions, condensation can occur in the air tube. If this occurs, after the treatment operate the unit with the tube connected but without the nebulizer to dry the tube.

5. Assemble clean nebulizer and add prescribed medication.
6. Attach mouthpiece or mask onto medication cup output port.

NOTE- please use in conjunction with compatible nebulizer kit originally supplied by Trimpeks combination with other nebulizer kit can affect flow characteristics and degrade system performance.

7. Connect the other end of the air tube to the nebulizer air-inlet connector.
8. Press the power switch “On”(I) to start the compressor.
9. Begin treatment by placing the mouthpiece between lips. With mouth closed, inhale deeply and slowly through mouth as aerosol begins to flow, hold for 2-3 seconds, then exhale slowly through the mouthpiece. If treatment needs to be interrupted, simply press power switch “Off” (O).

NOTE- your physician may recommend a specific breathing pattern during treatment. If so, follow his/her recommendations.







10. If an aerosol mask is used, place the mask over mouth and nose. As aerosol begins to flow, inhale deeply and slowly through mouth, then exhale slowly.

DANGER

Do not expose eyes to aerosol. Inhaled medications delivered through mask may inadvertently deposit in the eyes and result in eye irritation.

11. When treatment is completed, turn unit off by pressing the power switch to the “Off”(O) position, unplug the unit from the wall power outlet.

NOTE- stop nebulizer therapy at, or one minute after, the onset of sputter.

NEBULIZER CLEANING

All parts of nebulizer, except air tube, should be cleaned according to the following procedure. Your physician may specify other specific cleaning procedure. If so, follow their recommendations.

WARNING

To prevent possible risk of infection from contaminated medication, cleaning of the nebulizer is recommended after each use. Disinfecting is recommended once per day.

Clean After Each Use:

1. With power switch in the "Off"(O) position, unplug power cord from wall outlet.
2. Disconnect air tube from the air-inlet connector and set aside.
3. Disassemble mouthpiece or mask from medication cup. Open nebulizer by turning baffle counter-clockwise and removing.
4. Wash all items, except air tube, in a warm water/dishwashing detergent solution. Rinse under warm tap water for 30 seconds to remove detergent residue. Allow to air dry.

Disinfect daily:

1. Using a clean container or bowl, soak items in three parts warm water (55-65°C) to one part white vinegar (>5% concentration) for 60 minutes.
2. With clean hands, remove items from disinfectant solution, rinse under warm tap water, and air dry on a clean paper towel.

NOTE- *do not towel dry nebulizer parts, this could cause contamination.*

WARNING

To prevent possible risk of infection from contaminated cleaning solutions, always prepared fresh solution for each cleaning cycle and discard solution after each use.

Keep the outer surface of the air tube dust-free by wiping regularly. Air tube does not have to be washed because only filtered air passes through it.

NOTE- *the reusable nebulizer kit will last up to 6 months. Proper cleaning will also help extend the life of the nebulizer kit.*

NOTE- *ageing of material of nebulizer kit may occur after number of cleaning and/or disinfection cycles. If you find any discolour or physical change, please replace nebulizer kit with a new one.*

COMPRESSOR CLEANING

1. With power switch in the "Off" (O) position, unplug power cord from wall outlet.
2. Wipe outside of the compressor enclosure with a clean damp cloth every few days to keep dust-free.

DANGER

Do not submerge in water; doing so will result in compressor damage.

FILTER REPLACEMENT

1. Filter should be replaced every 30 days or sooner if filter discolours. This depends on the quality of ambient air.
2. Remove filter by grasping it firmly and pulling out of the unit. Discard the old filter.
3. Replace with a new one. Additional filters should be purchased from your dealer.
4. Push new filter into position.

CAUTION- *do not reuse or wash filter or substitute any other material such as cotton for an air-inlet filter. It will result in compressor damage.*

MAINTENANCE

All maintenance must be performed by a qualified provider or authorized service centre.


DANGER



Electrical shock hazard. Do not remove compressor enclosure. All disassembly and maintenance must be done by qualified provider.

SPECIFICATIONS

Power	AC 230V/50Hz
Power consumption	0,7 A
Max. fill volume	6 ml max.
Sound level	≅ 65 dB at 1 meter
Droplet distribution*	Refer to charts below (based on light-scattering technique)
Nebulization rate*	≅ 0.2ml/min
Residual volume	≅ 1,0ml
Max. compressor pressure	200 to 300 kPa
Operating pressure*	70 to 106kpa
Compressor air flow*	6 ~ 8 lpm
Operating temperature	5°C to 40°C (50°F to 104°F)
Operating Humidity	up to 95% RH
Storage/Transport Temperature	-20°C to 70°C (-4°F to 158°F)
Storage/Transport Humidity	10 to 95 % RH
Dimension (L x W x H)	294 x 170 x 100mm (9.7" x 5.6" x 3.3")
Weight	4.3lbs (2.0kg)
Standard Accessories	Medication cup, Air Tube, Mouth piece, Filters (5pcs), Adult mask, Child mask
Optional Accessories	Carrying bag
Service life	6 months (nebulizer kit); 1,400 hours (compressor)

*Conditions may vary based on the altitude, barometric pressure and temperature.

Equipment Classifications:

Protection against electric shock:	Class II
Degree of protection against electric shock:	Type BF
Mode of operation:	Continuous operation

Compressor is not suitable for use in the presence of flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

DISPOSAL

Please dispose of the device in accordance with EU Directive 2002/96/EC - WEEE.
 If you have any question, please contact the local authorities responsible for waste disposal.

GUIDELINE AND DECLARATION

WARNING

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical electrical equipment.

The equipment should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.

NOTE- The EMC tables and other guidelines provided information to the customer or user that is essential in determining the suitability of the Equipment for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment to perform its intended use without disturbing other Equipment or non-medical electrical equipment.

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The Compressor Nebulizers are intended for use in the electromagnetic environment specified below. The user of the devices should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Compressor Nebulizers use RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interface in nearby electronic equipment.
RF emissions CISPR 11	Class B	The devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliant	

Guidance and manufacturer's declaration – Electromagnetic Immunity			
The Compressor Nebulizers are intended for use in the electromagnetic environment specified below. The user of the Compressor Nebulizers should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact $\pm 2, 4, 8, 15$ kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5, 1$ kV line(s) to line(s)	± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions	$<5\% U_T$ ($>95\%$ dip in U_T)	($>95\%$ dip in U_T) for 0.5 cycle	Mains power quality should be that of a typical commercial or

and voltage variations on power supply input lines IEC 61000-4-11	for 0.5 cycle 70% U _T (30% dip in U _T) for 25cycles <5% U _T (>95% interruption in U _T)	(>95% dip in U_T) for 1 cycle (30% dip in U_T) for 25cycles (>95% interruption in U_T) for 250cycles	hospital environment. If the user of the Compressor Nebulizers requires continued operation during power mains interruptions, it is recommended that Compressor Nebulizers be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	If image distortion occurs, it may be necessary to position Compressor Nebulizers further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficient low.
NOTE U _T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – Electromagnetic Immunity			
The Compressor Nebulizers is intended for use in the electromagnetic environment specified below. The user of the Compressor Nebulizers should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V_{rms}	Portable and mobile RF communications equipment should be used no closer to any part of Compressor Nebulizers , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{V1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3, 10 V/m 80 MHz to 2700MHz	10 V/m	$d = \left[\frac{3,5}{E1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E1} \right] \sqrt{P}$ 800 MHz to 2700MHz Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

		<p>and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations, electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>		
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM And FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Compressor Nebulizers is used exceeds the applicable RF compliance level above, the Compressor Nebulizers should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Compressor Nebulizers.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m</p>		


Recommended separation distances between portable and mobile RF communications equipment and the Compressor Nebulizers			
<p>The Compressor Nebulizers are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Compressor Nebulizers can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Compressor Nebulizers as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2700MHz
	$d = \left[\frac{3,5}{V1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E1} \right] \sqrt{P}$	$d = \left[\frac{7}{E1} \right] \sqrt{P}$
0.01	0,12	0,12	0,23
0.1	0,38	0,38	0,73

1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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