plusmed

health management

Model: **pM-KN01**

OKSİJEN KONSANTRATÖRÜ
OXYGEN CONCENTRATOR

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EN Instruction Manual

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Kalite garantisi

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EN Foreword

Thank you for purchasing our products, hoping you will be satisfied with our products.

This operation manual contains function, operation steps, attention, basic trouble solution and so on.

To ensure your efficient use of the machine, please have a close read of this operation manual before operating it.

Maybe, there are some pictures, which are different from what you have seen in the real model in this manual.

Safety notice

This product cannot be used for life saving, it is suggested that if any patient who needs oxygen treatment, please follow doctor's advice to choose the right flow and period for oxygen before using the machine.

Please note the following special statements, used throughout this manual, and their significance:

Note: Explanatory information.

 Damage: Action could result in damage to equipment.

 ATTENTION: Action could result in personal injury.

 WARNING: Action could result in fire and explosion.

 ATTENTION: Consult Instructions for Use.

 Electrical category: CLASS II

 Electrical category: TYPE B

 : The date of Manufacturer

 : The Manufacturer

 0197: Indicating its conformity with the Medical Devices Directive 93/42/EEC. The 1984 is the number of the Notified Body.

 No open flame: Fire, open ignition source and smoking prohibited

 Do not Smoke

 Separate collection for electrical and electronic equipment.

Before Installation

 Damage: Only use stable and safe electrical power sources.

 Damage: If the electrical power source becomes unstable, discontinue use.

Placement

 **ATTENTION:** The oxygen concentrator should be set to use in an environment without dust, corruption or toxicological harm gas.

 **ATTENTION:** The oxygen concentrator should be located in well-ventilated space, in case there are polluted air or smog in the oxygen.

 **ATTENTION:** Do not place items on top of the concentrator.

 **ATTENTION:** Always place the concentrator on a hard surface. Never place the concentrator on a surface such as bed or couch, where the concentrator may tip or fall.

Fire Warning

 **WARNING:** For oxygen can be combustion-supporting, keep oxygen concentrator far away from naked light or fire resource, no smoking or naked light around the patient.

 **WARNING:** Keep the concentrator away from flammable and explosive areas.

 **WARNING:** Textiles and other materials that normally would not easily burn, can ignite and burn with great intensity in oxygen enriched air.

 **WARNING:** A spontaneous and violent ignition may occur if oil, grease or greasy substances come in contact with oxygen under pressure. ALWAYS keep these substances away from the oxygen concentrator.

Maintenance

 **ATTENTION:** Carrying out maintenance the device for every 2500 hours is recommended

 **ATTENTION:** Before cleaning the dust on the net of oxygen concentrator, the plug must be pulled out in case electric attack

 **ATTENTION:** The humidifier, filter cotton and filter are the items needed to clean, among which humidifier should be cleaned every 3 days, and external filter cotton (coarse dust filter) should be cleaned every 100 hours, and internal filter (air filter) should be changed every 500 hours.

Radio Frequency Interference

Most electronic equipment is influenced by Radio Frequency Interference (RFI). When there is strong electromagnetic interference, maybe the LCD will be slightly affected, but the machine is still running. ALWAYS exercise CAUTION with regard to the use of portable communications equipment in the area around such equipment.

Additional Safety Warnings

 **WARNING:** Do not put nasal tube under bed or cushion, the oxygen that caused by machine turning on without breathing may be combustion-supporting.

 **WARNING:** DO not reach for a concentrator that has fallen into water. UNPLUG IT IMMEDIATELY.

 **WARNING:** NEVER leave the concentrator unattended when plugged in.

 **ATTENTION:** ALWAYS supervise closely when this product is used near children or those who require close super-

vision.

⚠ ATTENTION: If any adverse reaction appeared or revealed during taking oxygen , please contact with equipment supplier or doctor as soon as possible.

⚠ ATTENTION: For serious patients, set an indicating device additional, any adverse reaction appeared, please contact with equipment supplier or doctor as soon as possible.

⚠ ATTENTION: Turn off the switch if nobody takes oxygen.

⚠ ATTENTION: In using the machine, do not open the front and back cover at all. In case there are quality problem, do not dismantle it secretly. Any alarm or other abnormal phenomenon has been found, connect with equipment supplier or factor.

⚠ ATTENTION: Ensure the bottom smooth exhaustion during operating, or else the machine will be over-heated.

⚠ ATTENTION: There is intermitted exhaustion sound during operating (13 seconds in intermission).

⚠ ATTENTION: 5 minutes are needed from oxygen concentrator from warming up to reach regular function.

⚠ ATTENTION: The machine is only for oxygen supply, and the oxygen concentration will be up to 90% when air outlet reaches its nominal flow.

⚠ ATTENTION: Humidifier shall adopt distilled water or cold boiled water, added water shall be kept under the scale line.

⚠ ATTENTION: Use the humidifier together with the machine, do not replace it at will, or else it may cause patient uncomfortable or other harms.

⚠ ATTENTION: In case the indicator shows abnormal oxygen, operator should declare to dealer or factory in favor of maintain.

⚠ ATTENTION: Once open the adjust knob for flow in full, but flow meter shows zero, turn off the machine immediately and have a check for trouble.

⚠ ATTENTION: Do not turn on or off frequently, To restart the machine after turning off, no less than 5 minutes are necessary (namely, exhaust internal gas of the machine completely, for if air compressor turns on with pressure, its life will be shortened)

⚠ ATTENTION: Turn up the flowmeter switch immediately when power switch is turned on.

⚠ ATTENTION: Refresh the water in the dampen bottle every 2-3 days, especially in summers. If do not use it in several days, please pour out the water completely, and wipe dry the bottle.

⚠ ATTENTION: Use the oxygen tube and humidifier together with the machine or those of the same model, if change to use other model devices, please ensure close connection with the oxygen generator. The absorbing oxygen tube is only for the patient, and do not junk it at will.

⚠ ATTENTION: The oxygen tube, oxygen mask that have touched with the patient should keep clean, disinfected and sterilized.

⚠ ATTENTION: The oxygen tube that have touched with the patient after each operation should be disinfected by wiping it with 75% medical use alcohol or other disinfecting methods. To prevent cross infection, do not share oxygen tube.

Safety notice

 ATTENTION: This unit is not a life-support device, and in certain circumstances oxygen therapy can be hazardous, it is suggested that if any patient who needs oxygen treatment, please follow doctor's advice to choose the right flow and period for oxygen before using the oxygen concentrator.

 ATTENTION: In the event of an alarm, you observe your oxygen concentrator not working properly, or if you feel discomfort, consult your Equipment Provider and /or your physician immediately.

 ATTENTION: Use only voltage specified on rating label.

 ATTENTION: This device manufactures high concentration oxygen, which promotes rapid burning.

Keep oxygen concentrator far away from open flames and no smoking around the patient.

 ATTENTION: Do not leave a nasal oxygen cannula under bed coverings or chair cushions. If the unit is turned on without use, the oxygen will help the flammable material get fire.

 ATTENTION: Use no lubricants, grease, or petroleum-based products on or near your oxygen concentrator.

 ATTENTION: Electrical shock hazard. Do not remove covers while the unit is plugged in. Only your Equipment Provider or a qualified service technician should remove the covers or service the unit.

 ATTENTION: Care should be taken to prevent the unit from getting wet or allowing water to enter the unit.

 ATTENTION: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

 ATTENTION: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the JAY-5 medical oxygen concentrator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

 ATTENTION: The oxygen concentrator should be set to use in an environment without dust, corruption or toxicological harm gas.

 ATTENTION: Do not place the oxygen concentrator in surroundings where its airflow is obstructed.

 ATTENTION: Do not place items on top of the concentrator.

 ATTENTION: Always place the concentrator on a hard surface. Never place the concentrator on a surface such as bed or couch, where the concentrator may tip or fall.

 ATTENTION: NEVER leave the concentrator unattended when plugged in.

 ATTENTION: Ensure the bottom smooth exhaustion during operating, or else the oxygen concentrator will be overheated.

 ATTENTION: 5 minutes are needed from oxygen concentrator from warming up to reach regular function and nominal performance.

NOTE: If oxygen does not seem to flow, first verify that the flowmeter ball is registering a flow. Then, place the tip of the cannula into a glass of water; if bubbles come out of the cannula, oxygen is flowing. If bubbles do not appear, turn off the oxygen concentrator immediately and refer to Troubleshooting.

NOTE: There is never a danger of depleting the oxygen in a room when you use your oxygen concentrator.

Requirement of environment protection

The materials used in the system won't create environment hazard. The packing materials of the system are recyclable, and they must be collected and disposed according to the related regulation in the country or region where the package of the system or its accessories is opened. The nasal oxygen tube is made of medical PVC, and if it is thrown away, it could not be bio-degradable, so it will cause the pollution. Any material of the system that may cause pollution in the environment, must be collected disposed strictly complied with the local rules and requirements.

Product introduction

pM Series oxygen concentrator model adopt pressure swing adsorption principle, which can separate oxygen, nitrogen and other gas from the air, at constant temperature, as soon as power is connected, the oxygen that meets medical use standards can be separated from air constantly. Oxygen is generated by pure physical method. The generator can supply ONE patient simultaneously, with adjustable and stable oxygen flow, safe and reliable, low cost. The key parts of the generator adopt anti-tiring and anti-aging design, and the planned life of the whole generator reaches up to 20,000 hours. There is no influence on indoor oxygen percent during the generator operating.

Using condition

Ambient temperature: 10°C-40°C

Relative humidity: 30%-85%

Air pressure: 700 hPa-1060 hPa

No corrosive gas and strong magnetic field around.

Altitude: Up to 2286m without degradation; Consult your equipment provider for further information regarding to 2286m to 4000m

Scope of application:

1. For Medical Use

oxygen supplied by the concentrator is beneficial to cure the disease or heart and blood vessel system, chronic pulmonary system, the brain and blood vessel system, chronic pulmonary tuberculosis, and other oxygen lacking symptoms, etc.

2. For Health care

oxygen can be used for athletics and intellectuals and brainworkers, etc. to eliminate fatigue and also suit for the departments of health care, sanatorium, healthy, plateau military camps and hotels and other places where need oxygen.

3. Device can be used by pregnant, nursing women and children under 14 years old.

Accessories:

Your concentrator includes the following components:

- Intake air filter (two pieces, part number:GL-01)
- Secondary filter (one piece, part number:GL-02)

The concentrator comes with two air filters and one secondary filter already installed.

Technical Parameter

1) Model description: pM-KN01

Model	pM-KN01
Rated power (VA)	550
Operation voltage (V/Hz)	AC230±10% 50±1
Oxygen flow (L/mln)	0-5
Oxygen concentration (%)	93%±3%
Outlet pressure (Mpa)	0.04—0.07
Alarm	Power failure low & high pressure ;
Sound level (dBA)	≤50
LCD display	accumulating timing; present timing; timing;
Large LCD display	Switch times; pressure digital(accuracy:0.001MPa); accumulating timing(range:0-10000hours); present timing(accuracy:1 minute); presetting timing(accuracy:1 minute)
Electrical category:	Class II Type B
Product category:	Class II a
Net Weight (Kg)	26
Fuse	T6.3AL/250V

2) Nasal oxygen tube

Nasal oxygen tube is made up of intubation, oxygen tube, inlet oxygen tube, movable gripper tube and nasal tube. And it is made of medical PVC. There are two tips, over ear style. It must be flexible hose, transparent and clean. There is no impurities and no kink. When inputting the 50Kpa pressure the nasal oxygen tube must be non-disclosure continued 15S. The various components can withstand vertical stress of 20N in 15S without break. With Ethylene oxide sterilization, since the date of sterilization, the nasal oxygen tube could be saved for three years.

3) Humidifier

The humidifier must be colorless and transparent or translucent. And it is marked with the maximum and minimum liquid level instruction line. It should not be broken when it withstands the pressure of not less than 0.4Mpa.

Name and Function

Diagram 1

Diagram 2

1) Indicating Lamp

As showing in Diagram 1:

Total 6 indicating lights, and the code of them are 1,2,3,...6 from left to right, and their indication for each model are as follows:

P.O.: power switch (green lamp)

P.F.: power failure (red lamp)

L.P.: low pressure (yellow lamp)

H.P. (T): high pressure (red lamp) / overheated temperature (red lamp)

L.O2: oxygen purity is < 85 %, (red lamp)

H.O2: oxygen purity is \geq 85%, (blue lamp)

2) Power switch

3) Oxygen flow meter

The location of float in the oxygen flow meter shows the outlet oxygen flow (L/min.).

4) Knob of oxygen flow meter switch

The other name of knob of oxygen flow meter switch is flow control valve. It adjusts and controls the outlet oxygen flow. Do not rotate it over-forced, or else it is easy to damage the valve core. Rotate it counterclockwise to turn on, clockwise to turn off.

5) Intake air filter

Once need to replace, replace it by special filter for the machine.

6) Back label

7) LCD display (Liquid crystal display)

a.It can display the following diagram 3.

b.At the beginning of starting the machine, the LCD screen is bright, after 15 minutes it will return to screen saver mode. But if you press the right key during working, the screen will bright again.


Diagram 3

8) Timing keys

As showing in Diagram 1: “+”.”-” These two keys used for timing. Press “+”button, increase timing by 10 minutes step. Press “-”button, decrease timing by 10 minutes step. When press“-” button till 0 minute, machine will turn off by itself, at this time press the “+”button, machine will restart.

9) Humidifier

Humidifier which used for humidifying oxygen and preventing throat and nasal mucosa stimulated by dry oxygen and dry hard sputum difficult to spit out.

10) Power socket

11) Fuse

12) Outlet for Atomization

13) Knob of Atomizer switch

Operation steps

1. As showing in Diagram 4: Pull off the silicone tube as black arrows direction ,which is connected with oxygen output ;then take down the humidifier bottle along clockwise direction of the red arrow, Fill in with distilled water or cooling boiled water Between min scale and max scale, then screw down the bottle. At last put the silicone tube back connect to the oxygen output.

Diagram 4

2. Connect the power, put the plug of power line connected with the power socket of the oxygen generator, and the other end of the plug connects with indoor power socket, turn on the power switch. (See diagram 5)


Diagram 5

3. Adjust well oxygen output flow according to the request, the red scale of the flow meter is MAXIMUM recommended oxygen absorbing flow (counterclockwise—on, clockwise—off, for common single treatment: pM-KN01 adjusted to 5L/min; For double use: pM-KN01 adjusted to 3-5L/min.(see diagram 6)


Diagram 6

4. Insert the intake end of absorbing oxygen tube onto the outlet of humidifier, then set the absorbing oxygen tube over patient's ears, insert the nasal tube into patient's nostrils to absorb oxygen; the best absorbing time for health care keeps 40-50 minutes per time, absorbing time for medical treatment shall be followed doctor's advice. (see diagram 7)


Diagram 7

 **WARNING:** 5 minutes are needed from oxygen concentrator from warming up to reach regular function and nominal performance.

5. After breathing oxygen, turn off the power, and pull up the nasal oxygen cannula from humidification bottle. If turn off the machine without pull off the nasal oxygen cannula, it may cause oxygen tube bending and made remain water in the humidification bottle backflow to machine, further lead to machine failure working .When discontinuous using, please unplug the power plug.(see diagram 8)

Atomization operation methods

The atomization function of the concentrator is applicable to help to cure the sufferers of chronic pulmonary tuberculosis and respiratory system, etc.

⚠ :Atomization treatment must be under advice and suggestion of your physician.

⚠ :Using distilled water to do the atomization for several seconds after each operation may lighten the crystallization caused by medical solution.

⚠ :If atomizing cannot work then open the cover of the bottle and add clean water in small amount. Rotate the white ball which lies in the bottle with the gas resource connected and select the proper angle to gain a better atomization.

- a. Open the cover of medical cup, and add atomization remedy that needed, then close the cover.
- b. Connect the joint of atomization nozzle (or mask) with the cover of medical cup, and then connect the other end of atomizer connection tube with the atomization outlet and turn on atomizer.
- c. Turn on the power of oxygen concentrator, and shut up flow meter, then it is ready for atomization treatment.
- d. Do clean the atomization devices followed by the instructions of the atomization devices.

Unplug the power plug


Diagram 8

6. If the patient need timing oxygen absorbing, please refer to instruction on page 7 item 9 (Timing keys).

Maintenance

1. In the condition of power off, make a clean for the outside body by soft towel with little detergent, and then wipe it up with dry towel, once or twice per month.(see diagram 9)

Diagram 9

2. It is a critical step for daily maintenance to clean intake air filter, at least twice a month.

Detail steps: take off the two intake air filter on both sides of the body, clean them with detergent and clean it out with clean water completely, get ride of extra water and dry up naturally, finally set back after dry up.(see diagram 10)

air filter

Diagram 10

▲ Damage: Do not operate the concentrator without the filters installed, or while filters are wet. These actions could permanently damage the concentrator.

3. Change of Secondary Filter (Under the humidifier)

First of all remove the intake filter and then turn off the secondary filter counterclockwise as shown in Figure 11. Remove the filter cloth and change it with the new filter.


Diagram 11.1

Diagram 11.2

▲ Damage: Do not operate the concentrator without the filters installed, or while filters are wet. These actions could permanently damage the concentrator.

4. Replacement of fuse tube

Take off the cover of fuse, which is on the power socket, dismantle the fuse tube off by small screwdriver. Close the cover of fuse after fuse tube is replaced.

The other fuse tube is located at the intake of internal power line; the method of replacement is the same with that above. (see diagram 12)



Diagram 12

5. The battery of pM Series oxygen concentrator model is NI-MH charge battery (3.6V/40mAh), which welded on the main board and there is a charge circuitry in the board, so it could be used by a long time. But if you could not use the oxygen concentrator continually, the battery will be at the dormancy state and its life will be shorter. So, please charge the battery one time every month if you could not use the oxygen concentrator continually.

6. If the user want the circuit diagram the list of critical components the detail of maintenance and repair, we will supply the data about all the repairable parts which we can supply.

Troubles and solution

No.	Trouble	Causes	Solution
1	No operation after power connected	<ol style="list-style-type: none"> 1. No connection between circuit of oxygen generator and power 2. Circuit of fuse protector broken 3. Capacitor of compressor broken 4. Compressor broken 	<ol style="list-style-type: none"> 1. Check out whether switch, plug, power line in good connection. 2. Replace the fuse protector and find the cause 3. Replace start capacitor 4. Have the compressor replaced
2	No oxygen out or tiny outtake flow	<ol style="list-style-type: none"> 1. Folded inside oxygen tube, no smooth outtake 2. Filter clogged, no smooth intake 3. The cover of dampen bottle leaking 	<ol style="list-style-type: none"> 1. Connect the oxygen tube again 2. Clean the filter 3. Take off the cover, screw well the cover, block the outtake by thumb after turning on, and there will some sound from the humidifier after 5 second around (the safety valve of humidifier turns on)

3	No exhaust sound	1. Air controller cannot work 2. Electrical control board cannot work	1. Have air control valve replaced 2. Have electric control board replaced
4	Too noisy exhaust	1. The joint of exhaustion muffler fallen off 2. Exhaustion muffler broken	1. Connect the joint well 2. Have the muffler replaced

EMC Declaration

The pM-KN01 needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents;

Portable and mobile RF communications equipment can affect the pM-KN01.

All cables and maximum length of cables, Transducers and other accessories with which the manufacturer of the pM-KN01 claims compliance with the requirements, Accessories that do not affect compliance with the requirements of these sub clauses need not be listed. Accessories, transducers and cables may be specified either generically or specifically.

NOTE

Transducers and cables sold by the manufacturer of the pM-KN01 as replacement parts for internal components need not be listed.

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of The pM-KN01 as replacement parts for internal components, may result in increased emissions or decreased immunity of The pM-KN01.

Guidance and manufacturer's declaration – electromagnetic emissions		
The pM-KN01 is intended for use in the electromagnetic environment specified below. The customer or the user of the pM-KN01 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The pM-KN01 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The pM-KN01 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This pM-KN01 is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the pM-KN01 or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The pM-KN01 is intended for use in the electromagnetic environment specified below. The customer or the user of the pM-KN01 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance

Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output Lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the pM-KN01 requires continued operation during power mains interruptions, it is recommended that the pM-KN01 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable Note: The pM-KN01 does not contain components susceptible to magnetic fields, such as Hall elements or magnetic field sensors. Therefore, the EUT is deemed to meet the requirement without actual testing.	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the pM-KN01, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.1 \sqrt{P}$ $d = 1.1 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the pM-KN01 is used exceeds the applicable RF compliance level above, the pM-KN01 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the pM-KN01.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the pM-KN01

The pM-KN01 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the pM-KN01 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the pM-KN01 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2,5 GHz $d = 1.17 \sqrt{P}$
0.01	0.12	0.12	0.07
0.1	0.37	0.37	0.22
1	1.17	1.17	0.70
10	3.69	3.69	2.21
100	11.67	11.67	7.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Condition for transportation and storage

Environment temperature scale: -20°C-45°C

Comparative humidity scale: ≤95%

Air pressure scale: 500 –1060 hPa

Quality Warrant

The warranty period for the product is 2 years.

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