




## UltraScan II

INFRARED THERMOMETER  
EAR & FOREHEAD SCAN

- 1 Measured on forehead or ear
- 2 Fast reading
- 3 10 reading memory
- 4 3-color Backlight
- 5 Result in 1 second
- 6 Automatic power-off



Operating Condition	Operating Temperature : 10°C~40°C(50~104°F) Relative Humidity : 15%~85% Atmospheric Pressure : 70-106Kpa
Transport and Storage Condition	Storage Temperature : -25~55°C / -13~131°F Relative Humidity : 15%~95% Atmospheric Pressure : 70-106Kpa
Operating Mode	Adjusted Mode (Body Mode)